

QuantiGene[®] Plex Assays for DNA Copy Number

Get highly sensitive and accurate multiplex copy number information for any DNA sample

Ultimate precision and accuracy with QuantiGene[®] Plex DNA Assays

- No DNA extraction or purification required
- No amplification
- Discriminate 0.5 to 4 copies
- Multiplex 3-34 targets
- Only 10 ng genomic DNA required

Robust across various sample types

- Formalin-fixed, paraffin embedded (FFPE) samples
- Fresh/frozen tissue
- Blood
- Cell lysates

Use QuantiGene DNA assays to confirm copy number variations (CNVs)

- CNVs from microarrays
- Cancer biomarkers
- Chromosomal breakpoints
- Areas of amplification/deletion
- Aneuploidy

Choose from off-the-shelf QuantiGene[®] Assays or customize your own in three to four weeks.

